



Biomark

PRE-LOADED TRAYS

HIGH PERFORMANCE PIT TAGS



Biomark **Pre-loaded Trays** has tags loaded in single use needles and packaged in a tray of 100. The tags are secured inside the needle, so no worrying about tags falling out before being implanted. It also eliminates the need to disinfect between implants, making it a very efficient and safe tagging system.

PIT TAG IMPLANTERS

Compatible with the MK10 syringe with N125 needle, the MK7 Combo, Biomark Pre-load System and Biomark Pre-load Sterile



FEATURES

- Enhanced read range performance
- Low-frequency 134.2 kHz operation
- 64-bit identification code
- ISO 11784/11785 FDX-B compliant
- Biocompatible glass encapsulation


APPLICATIONS

- In-vivo Animal Identification
 - Fisheries (marine & fresh water)
 - Small & Large Mammal
 - Reptiles & Amphibians
 - Birds & Bats
- Generic Object Identification
 - Rocks
 - Trees
 - Plants

GPT12™ Pre-load Tray
 HPT9™ Pre-load Tray
 HDX12™ Pre-load Tray



APT12™ Pre-load Tray



MiniHPT8™ Pre-load Tray
 MiniHPT10™ Pre-load Tray

ALL BIOMARK PIT TAGS TRAYS ARE 134.2 KHZ, ISO 11784/11785 COMPLIANT AND ICAR APPROVED.

Tag Model	Bulk	Pre-load Trays	Individual Pre-load	Sterile Syringe/ Needle Combo	Compatible Needle/Syringe
MiniHPT8 (8.4 mm x 1.4 mm)	✓	✓	✓	✓	N165 MK165 / MK65
HPT9 (9 mm x 2.1 mm)	✓	✓	✓		N125 MK10/MK7 / MK25
MiniHPT10 (10.3 mm x 1.4 mm)	✓	✓	✓	✓	N165 MK165 / MK65
APT12 (12.5 mm x 2.03 mm)	✓	✓	✓		N125 MK10/MK7 / MK25
GPT12 (12.5 mm x 2.1 mm)	✓	✓	✓	✓	N125 MK10/MK7 / MK25
HDX12 (12 mm x 2.12 mm)	✓	✓	✓		N206 MK10



Biomark 

Scan the code to access this and other resources in the Biomark Digital Library, or visit: biomark.com/library

 **MERCK**
Animal Health

© 2021 Intervet Inc., d/b/a Merck Animal Health, a subsidiary of Merck & Co. Inc., All rights reserved.